N⊗rth **BULK CARGO LIQUEFACTION** THE CAN TEST POCKET CARD

For more information, please visit www.nepia.com 😏 😭 in 🖸 Copyright © 2019 The North of England P&I Association Limited

THE CAN TEST



 A rigid cylindrical can - (0.5 to 1 L capacity).
Half fill with a sample of the material.
Bring it down sharply from a height of about 0.2 m
Repeat 25 times - at 1 - 2 second intervals.
Take high definition pictures and video
Examine the surface for free moisture or fluid conditions.
If free moisture or a fluid condition appears DO NOT LOAD and SEEK ADVICE.
www.nepia.com/bitesize-training

For further information please contact our Loss Prevention team on loss.prevention@nepia.com